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
**USER'S GUIDE**

**MANUAL DE INSTRUCCIONES**

**MODE D'EMPLOI**

Abnil 2012

PT | INDICE

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
**While working under the fume cupboard, the following safety instructions should be observed:**

- Immediately open all the windows whenever it smells like gas. No electric switches should be turned on.
- Immediately call the fire brigade in case of fire. Until the arrival of the fire brigade, the fire should be fought with the available fire extinguishers.
- Always keep the emergency exits unobstructed, as well as the access ways to the showers and the eye-washers.
- Always wear the appropriate protective equipment and apparel.
- In case of accident with electric power, immediately activate the EMERGENCY-OFF switch.
- Hazardous tasks should be performed by one person only, and in emergency cases another person should be available in order to offer first aid support.
- Any anomalies that might be detected should be immediately eliminated. The event should be reported to the senior member of staff.
- The work place can only be abandoned if a continuous control of ongoing tests is not required or in case any other person, with appropriate training, takes up the control of the tests.
- Any tests which cannot be interrupted by the end of normal working hours should only go on, without continuous control, if no accident can occur due to the safety measures already applied, or in case it is not possible to carry out the test on a different day.
- Hazardous substances should be stored in such a way as not to endanger human health or the environment.
- If the fume cupboard is used in a manner not specified by Laborial, the protection provided by the equipment may be impaired.

**Whenever working with explosive substances, the most restrictive standards should be applied.**

**2.1 - CONSTITUTION**

|    |                      |     |                      |     |                           |
|----|----------------------|-----|----------------------|-----|---------------------------|
| 1. | Upper body           | 6.  | Service handles      | 11. | Working area              |
| 2. | Lower body           | 7.  | Sliding window       | 12. | Upper light               |
| 3. | Service panel        | 8.  | Sliding glass panels | 13. | Deflecting/service panels |
| 4. | FC / FV01 Controller | 9.  | Window handle        | 14. | Upper front panel         |
| 5. | Power sockets        | 10. | Fixed window         | 15. | Collector                 |

### Sash window

The sash window is made by an aluminum frame, with two glass panels with transversal movement (sliding windows).

The work area becomes accessible when the window is opened. This window has an ascending movement for opening and a descending one for closing.

The horizontal sliding windows enable the operator to operate within the work area without having to fully open the window, thus physically protecting the body.

The maximum operational opening (500mm) is limited by a stoppage system built into the window.

The limiting break prevents the opening of the window to be higher than the maximum operational opening established by the standards. The break may be unlocked by an activation a device, thus enabling the window to open more than the operational maximum.

In case the hood is operating, an intermittent sound and visual alarm is activated whenever the opening of the window is greater than the maximum operational opening.

The sash window is equipped with an anti-fall system, in case any of the suspension cables fails.


A small collector in the internal part of the window handle drains and prevents any liquid spills from pouring out of the work area.

### Sash window Controller (Optional)

The sash window controller SC allows the window to close automatically, by resorting to an engine and a movement detector.

A photoelectrical barrier on the basis of the window detects any obstacles in the direction of the window, immediately interrupting its automatic closing.

### Command System

1. RS232 4 pole input - PC connection or for programming devices
2. Exhaustion ON/OFF Switch
3. Lighting ON/OFF Switch
4. Switch - Air flow in its minimum value (only applicable to the FV01 system)
5. LED Indicator - Sash window with opening above 500mm
6. Switch - Disconnect sound alarm
7. LED Indicator - Air flow below the minimum safety value
8. LED Indicator - Safe Air flow
9. Button - Air flow in its maximum value - (only applicable to the FV01 system)

### FC Control System

The FC control system allows monitoring the flow of air extraction through a pressure sensor.

Whenever the value of the air flow falls below the standard minimum safety value, a visual alarm and a continuous sound alarm are activated.

The sound alarm may be deactivated by pressing a button <sup>(6)</sup> on the control panel, while the visual alarm can only be deactivated after the problem is solved.

### FV01 Control System

The FV01 control system optional to the FC control system.

Additionally to the functionalities of the FC control system, this system enables control of the air flow (increasing or decreasing), manually or automatically.

This system makes it possible to compensate any fluctuations in the ventilation system.

Optionally, the extracted air flow can be controlled by the position of the window, in order to preserve a constant facial speed, in the section of the opening area. The air flow increases whenever the window opens, and decreases whenever the window closes.

Whenever the value for the air flow falls below the standard minimum safety value, a visual alarm and a continuous sound alarm are activated.

The sound alarm may be deactivated by pressing a button <sup>(6)</sup> on the control panel, while the visual alarm can only be deactivated after the problem is solved.

#### Lighting

The lighting system is designed in such a way that all the working area is illuminated with a lighting value above 450 lux.

Light projection is achieved through a luminary placed in the upper area, designed in such a way as to not dazzle the operator.

In case of any explosion, the luminary is released in order to relieve the pressure in the upper area, thus protecting the operator and the fume cupboard from any eventual damages.

Light control is made through a button on the control panels.

#### Service Panel

The service panel is developed so that it may contain the remote commands for service taps and valves, the power sockets and the air flow control panel.

#### Deflecting Panels

The fume cupboard is equipped with an upper horizontal deflective panel and with 4, 5 and 6 vertical deflective panels, for 1200, 1500 and 1800mm hoods respectively.



The vertical deflective panels shelter the taps and valves of the respective services to be installed in the hood (water, gases, etc).

#### Upper Front Panel

It is a removable panel in order to easily access the services installed in the upper area of the fume cupboard (collector, sash window sensor, controller, electric area).


## 2.2 - DIMENSIONS

- Height 2400mm
- Height with exhaust tube 2600mm
- Length 1200/1500/1800mm
- Depth 900mm
- Height of the work area 900mm
- Window maximum opening 500mm/650mm
- Useful working length 1158/1458/1758mm
- Useful working depth 675mm
- Useful working height 1275mm

## 2.3 – WORKING INSTRUCTIONS


### 2.3.1 – FC Control system

| Command/LED Indicator                      | Function                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 RS232 Input, 4 poles                     | Connection to a PC or to a programming device in order to program the control parameters.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| 2 ON/OFF Switch                            | <ul style="list-style-type: none"> <li>Only qualified technicians are authorized to change these parameters.</li> <li>Connects and disconnects the hood exhaust.</li> <li>Allows commissioning the hood controller. After the connection is done, the controller times about 5 seconds without monitoring any anomaly, in order to stabilize the ventilator.</li> <li>This button is green when connected.</li> <li>Connects and disconnects the light inside the Hood.</li> <li>Allows connecting and disconnecting the light inside the hood. This function is independent from the controller.</li> </ul> |
| 3 Lighting ON/OFF Switch                   | <ul style="list-style-type: none"> <li>Sets the exhaust flow in the pre-established minimum value.</li> <li>This button only works when the hood is ON.</li> <li>Issues a visual alarm whenever the window is positioned above 500mm.</li> <li>The visual alarm only fades out when the window is repositioned below 500mm.</li> </ul>                                                                                                                                                                                                                                                                       |
| (*) Switch – Air flow in its minimum value | Disconnects the sound alarms.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| LED Indicator – Window opened above 500mm  | <ul style="list-style-type: none"> <li>Allows disconnecting the sound alarms, regardless of the cause.</li> <li>Issues a visual alarm whenever the air flow is below the standard minimum safety value.</li> <li>The visual alarm only fades out when the air flow is greater than the minimum value.</li> </ul>                                                                                                                                                                                                                                                                                             |
| Sound alarm switch                         | Issues a luminous sign while the functioning is safe, i.e., whenever the air flow is above the standard minimum safety value.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| LED Indicator – Low air flow               | <ul style="list-style-type: none"> <li>Whenever the value is outside the range of correct parameters, this LED fades out.</li> <li>Sets the exhaust flow in the pre-established maximum value.</li> <li>This button only works when the hood is ON.</li> </ul>                                                                                                                                                                                                                                                                                                                                               |
| LED Indicator – Safe air flow              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| (*) Switch – Air flow in its maximum value |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

(\*) - Only applicable to the FV01 system

### 2.3.2 - FV01 Control system

| Command/LED Indicator                      | Function                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 RS232 Input, 4 poles                     | Connection to a PC or to a programming device in order to program the control parameters.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| 2 ON/OFF Switch                            | <ul style="list-style-type: none"> <li>Only qualified technicians are authorized to change these parameters.</li> <li>Connects and disconnects the hood exhaust.</li> <li>Allows commissioning the hood controller. After the connection is done, the controller times about 5 seconds without monitoring any anomaly, in order to stabilize the ventilator.</li> <li>This button is green when connected.</li> <li>Connects and disconnects the light inside the Hood.</li> <li>Allows connecting and disconnecting the light inside the hood. This function is independent from the controller.</li> </ul> |
| 3 Lighting ON/OFF Switch                   | <ul style="list-style-type: none"> <li>Sets the exhaust flow in the pre-established minimum value.</li> <li>This button only works when the hood is ON.</li> <li>Issues a visual alarm whenever the window is positioned above 500mm.</li> <li>The visual alarm only fades out when the window is repositioned below 500mm.</li> </ul>                                                                                                                                                                                                                                                                       |
| (*) Switch – Air flow in its minimum value | Disconnects the sound alarms.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| LED Indicator – Window opened above 500mm  | <ul style="list-style-type: none"> <li>Allows disconnecting the sound alarms, regardless of the cause.</li> <li>Issues a visual alarm whenever the air flow is below the standard minimum safety value.</li> <li>The visual alarm only fades out when the air flow is greater than such minimum value.</li> </ul>                                                                                                                                                                                                                                                                                            |
| Sound alarm switch                         | Issues a luminous sign while the functioning is safe, i.e., whenever the air flow is above the standard minimum safety value.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| LED Indicator – Low air flow               | <ul style="list-style-type: none"> <li>Whenever the value is outside the range of correct parameters, this LED fades out.</li> <li>Sets the exhaust flow in the pre-established maximum value.</li> <li>This button only works when the hood is ON.</li> </ul>                                                                                                                                                                                                                                                                                                                                               |
| LED Indicator – Correct air flow           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Button – Air flow in its maximum value     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

2.5 – Safety Rules for Operating the Fume Cupboard

- Never work with the fume cupboard, if the green LED Indicator is not ON.
- Keep the window closed during any test at all times.
- Open the transversal sliding window, only if absolutely necessary to intervene during the test.
- Try to keep the face as well as other parts of the body well protected by the transversal sliding window at all times during any intervention.
- Immediately close the window, in case you hear any acoustic sound and whenever a red alarm light goes on (in the control panel), during any test. Immediately inform everyone in the laboratory, as well as the senior staff. If possible, the test should be interrupted and the ventilation system checked.

3 – FLOW RATES

|                            | Fume Cupboard<br>1200 | Fume Cupboard<br>1500 | Fume Cupboard<br>1800 |
|----------------------------|-----------------------|-----------------------|-----------------------|
| Minimum flow value (Alarm) | 430 m3/h              | 550 m3/h              | 670 m3/h              |
| Project flow value         | 520 m3/h              | 660 m3/h              | 800 m3/h              |

- The data presented are minimum values and is therefore not recommended its use as basis values for the design of any ventilation system.

- The ventilation system connected to the fume cupboard have to be agreed with the Laborial SA

2.3.3 – Alarm visualization and recognition

| Event                                                             | Button/Light |        |       |       | Sound Alarm |
|-------------------------------------------------------------------|--------------|--------|-------|-------|-------------|
|                                                                   | Green        | Yellow | Red   | Green |             |
| Normal functioning                                                | ○            |        |       | ○     |             |
| low air flow                                                      | ○            |        | ○     |       |             |
| Window opening higher than the maximum operating opening (500mm). | ○            | ☆      |       |       | ☆           |
| Power failure – Controller connected                              |              |        | ☆ (1) |       | ○           |
| Power failure – Controller disconnected                           |              |        | ☆ (2) |       |             |

(○) – Continuous; (☆) – Intermittent; (☆) – Intermittent every 2 seconds

- The sound alarm may be deactivated by pressing the button 
- The sound alarm shall automatically disconnect after 10 minutes in case it is not deactivated manually.
- The visual alarm will only fade out after the issues have been solved:
  - Window opening below the safety limit (500mm);
  - Activation of power supply;
  - Safe air flow (above the standard minimum safety value).

2.4 – Forbidden Substances

It is forbidden to work with the following substances in the fume cupboard:

- Radioactive substances;
- Microorganisms;
- Genetically transformed substances.

4 – CLEANLINESS AND MAINTENANCE

- Regular cleanliness enhances safety and preserves the value of the equipment;
- All dirt should be immediately removed with the appropriate cleaning products and by using a soft damp piece of cloth;
- Any dirt which may be more difficult to remove should be eliminated with domestic cleaning products for commercial use, as long as they do not contain any abrasive or very alkaline components;
- Never use organic solvents for cleaning;
- Before starting any cleaning procedure, check whether there are any hazardous residual materials for your health in the work area.
- The controller is equipped with a rechargeable lead-acid battery 6V/1.2Ah. The battery must be exchanged at most once every 5 years. Contains substances harmful to the environment and should be collected and processed in the proper areas in accordance with CE standards in force

5 – TECHNICAL DATA

|                                     | 1200                         | 1500   | 1800   |
|-------------------------------------|------------------------------|--------|--------|
| <b>Exterior Dimensions</b>          |                              |        |        |
| Length                              | 1200mm                       | 1500mm | 1800mm |
| Depth                               | 900mm                        | 900    | 900mm  |
| Height                              | 2400mm                       | 2400mm | 2400mm |
| Height with exhaust tube (minimum)  | 2600mm                       | 2600mm | 2600mm |
| <b>Interior Dimensions (useful)</b> |                              |        |        |
| Length                              | 1158mm                       | 1458mm | 1758mm |
| Depth                               | 675mm                        | 675mm  | 675mm  |
| Height                              | 1275mm                       | 1275mm | 1275mm |
| <b>Working Height</b>               |                              |        |        |
| Window maximum opening (approx.)    | 900mm                        | 900mm  | 900mm  |
| weight (approx.)                    | 650mm                        | 650mm  | 650mm  |
| Operating Temperature               | 250Kg                        | 300Kg  | 350Kg  |
| Voltage                             | 0°C +55°C                    |        |        |
| Socket Protection                   | 230V/400V   50 HZ   Max. 30A |        |        |
|                                     | IP44                         |        |        |

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